

For Research Use Only

NBEAL1 Polyclonal antibody

Catalog Number: 25431-1-AP



Basic Information

Catalog Number:

25431-1-AP

Size:

150ul , Concentration: 750 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22125

GenBank Accession Number:

BC137214

GeneID (NCBI):

65065

UNIPROT ID:

Q6ZS30

Full Name:

neurobeachin-like 1

Calculated MW:

200 aa, 307 kDa

Observed MW:

280 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells,

Background Information

Neurobeachin-like 1 (NBEAL1) is a Golgi-associated protein that plays an essential role in cholesterol metabolism regulation.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 25431-1-AP (NBEAL1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.