

For Research Use Only

VSTM2L Polyclonal antibody

Catalog Number: 25457-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number: 25457-1-AP	GenBank Accession Number: BC033818	Purification Method: Antigen affinity purification
Size: 150ul, Concentration: 750 µg/ml by Nanodrop and 307 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 128434	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	Full Name: V-set and transmembrane domain containing 2 like	
Isotype: IgG	Calculated MW: 204 aa, 22 kDa	
Immunogen Catalog Number: AG22292		

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human liver tissue,
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yun Li	37030064	Biochem Biophys Res Commun	WB

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 25457-1-AP (VSTM2L antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 25457-1-AP (VSTM2L antibody) at dilution of 1:200 (under 40x lens).