

For Research Use Only

C1orf109 Polyclonal antibody, PBS Only

Catalog Number: 25552-1-PBS

Featured Product



Basic Information

Catalog Number:

25552-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22147

GenBank Accession Number:

BC018109

GeneID (NCBI):

54955

UNIPROT ID:

Q9NX04

Full Name:

chromosome 1 open reading frame
109

Calculated MW:

203 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25552-1-AP (C1orf109 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25552-1-PBS in a different storage buffer formulation.



WB result of C1orf109 antibody (25552-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-C1orf109 transfected LNCaP cells. This data was developed using the same antibody clone with 25552-1-PBS in a different storage buffer formulation.