

For Research Use Only

ZNF695 Polyclonal antibody, PBS Only

Catalog Number: 25556-1-PBS



Basic Information

Catalog Number:

25556-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22202

GenBank Accession Number:

BC023527

GeneID (NCBI):

57116

UNIPROT ID:

Q8IW36

Full Name:

zinc finger protein 695

Calculated MW:

515 aa, 60 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 25556-1-AP (ZNF695 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25556-1-PBS in a different storage buffer formulation.