

For Research Use Only

ARMC6 Polyclonal antibody

Catalog Number: 25558-1-AP



Basic Information

Catalog Number: 25558-1-AP	GenBank Accession Number: BC067125	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop;	GeneID (NCBI): 93436	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: armadillo repeat containing 6	
Isotype: IgG	Calculated MW: 501 aa, 54 kDa	
Immunogen Catalog Number: AG22278	Observed MW: 40-50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : COLO 320 cells,
Species Specificity: human	

Background Information

ARMC6 (Armadillo repeat-containing protein 6) is a Protein Coding gene. Diseases associated with ARMC6 include Leber Congenital Amaurosis 1. The function of this gene's protein product has not been determined.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25558-1-AP (ARMC6 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.