

For Research Use Only

DYDC2 Polyclonal antibody

Catalog Number: 25568-1-AP



Basic Information

Catalog Number:

25568-1-AP

Size:

150ul, Concentration: 750 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22250

GenBank Accession Number:

BC007374

GeneID (NCBI):

84332

UNIPROT ID:

Q96IM9

Full Name:

DPY30 domain containing 2

Calculated MW:

177 aa, 21 kDa

Observed MW:

21 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, pig

Positive Controls:

WB : pig oviduct tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



pig oviduct tissue were subjected to SDS PAGE followed by western blot with 25568-1-AP (DYDC2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.