

For Research Use Only

DOPEY1 Polyclonal antibody

Catalog Number: 25569-1-AP



Basic Information

Catalog Number:

25569-1-AP

Size:

150ul, Concentration: 800 ug/ml by Nanodrop and 360 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22244

GenBank Accession Number:

BC150619

GeneID (NCBI):

23033

UNIPROT ID:

Q5JWR5

Full Name:

dopey family member 1

Calculated MW:

2456 aa, 276 kDa

Observed MW:

277 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 25569-1-AP (DOPEY1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.