

For Research Use Only

DOPEY1 Polyclonal antibody

Catalog Number: 25569-1-AP



Basic Information

Catalog Number: 25569-1-AP	GenBank Accession Number: BC150619	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 800 µg/ml by Nanodrop and 360 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 23033	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	Full Name: dopey family member 1	
Isotype: IgG	Calculated MW: 2456 aa, 276 kDa	
Immunogen Catalog Number: AG22244	Observed MW: 277 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 25569-1-AP (DOPEY1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.