

For Research Use Only

DOPEY1 Polyclonal antibody, PBS Only

Catalog Number: 25569-1-PBS



Basic Information

Catalog Number:

25569-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22244

GenBank Accession Number:

BC150619

GeneID (NCBI):

23033

UNIPROT ID:

Q5JWR5

Full Name:

dopey family member 1

Calculated MW:

2456 aa, 276 kDa

Observed MW:

277 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS page followed by western blot with 25569-1-AP (DOPEY1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25569-1-PBS in a different storage buffer formulation.