

For Research Use Only

TMEM81 Polyclonal antibody

Catalog Number: 25634-1-AP



Basic Information

Catalog Number:

25634-1-AP

Size:

150ul, Concentration: 550 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22367

GenBank Accession Number:

BC061592

GeneID (NCBI):

388730

UNIPROT ID:

Q6P7N7

Full Name:

transmembrane protein 81

Calculated MW:

255 aa, 28 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, HepG2 cells, mouse brain tissue

Background Information

TMEM81 (also named HC3107 or KVLA2788) may be an integral component of the membrane.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25634-1-AP (TMEM81 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.