

For Research Use Only

TMEM81 Polyclonal antibody, PBS Only

Catalog Number: 25634-1-PBS



Basic Information

Catalog Number: 25634-1-PBS	GenBank Accession Number: BC061592	Purification Method: Antigen affinity purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 388730	
Source: Rabbit	UNIPROT ID: Q6P7N7	
Isotype: IgG	Full Name: transmembrane protein 81	
Immunogen Catalog Number: AG22367	Calculated MW: 255 aa, 28 kDa	
	Observed MW: 33 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

TMEM81 (also named HC3107 or KVLA2788) may be an integral component of the membrane.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25634-1-AP (TMEM81 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25634-1-PBS in a different storage buffer formulation.