

For Research Use Only

TARBP1 Polyclonal antibody

Catalog Number: 25647-1-AP



Basic Information

Catalog Number: 25647-1-AP	GenBank Accession Number: BC016023	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 860 µg/ml by Nanodrop;	GeneID (NCBI): 6894	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: TAR (HIV-1) RNA binding protein 1	
Isotype: IgG	Observed MW: 185 kDa	
Immunogen Catalog Number: AG22675		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25647-1-AP (TARBP1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.