

For Research Use Only

TPRG1L Polyclonal antibody

Catalog Number: 25663-1-AP



Basic Information

Catalog Number: 25663-1-AP	GenBank Accession Number: BC019034	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 µg/ml by Nanodrop;	GeneID (NCBI): 127262	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: tumor protein p63 regulated 1-like	
Isotype: IgG	Observed MW: 31 kDa	
Immunogen Catalog Number: AG22585		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25663-1-AP (TPRG1L antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.