

For Research Use Only

# CDKL1 Polyclonal antibody

Catalog Number: 25702-1-AP



## Basic Information

Catalog Number:

25702-1-AP

Size:

150ul , Concentration: 900 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22541

GenBank Accession Number:

BC104977

GeneID (NCBI):

8814

UNIPROT ID:

Q00532

Full Name:

cyclin-dependent kinase-like 1 (CDC2-related kinase)

Calculated MW:

358 aa, 42 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A549 cells, COLO 205 cells

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25702-1-AP (CDKL1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.