

For Research Use Only

C17orf64 Polyclonal antibody, PBS Only

Catalog Number: 25706-1-PBS



Basic Information

Catalog Number: 25706-1-PBS	GenBank Accession Number: BC048806	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 124773	
Source: Rabbit	UNIPROT ID: Q86WR6	
Isotype: IgG	Full Name: chromosome 17 open reading frame 64	
Immunogen Catalog Number: AG22769	Calculated MW: 236 aa, 27 kDa	
	Observed MW: 27-30 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

CHCT1 (CHD1 helical C-terminal domain-containing protein 1), also called C17orf64, may play a role in the regulation of apoptosis.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25706-1-AP (C17orf64 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25706-1-PBS in a different storage buffer formulation.