

For Research Use Only

ZNF419 Polyclonal antibody

Catalog Number: 25724-1-AP



Basic Information

Catalog Number: 25724-1-AP	GenBank Accession Number: BC008297	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 650 µg/ml by Nanodrop and 340 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 79744	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: zinc finger protein 419	
Isotype: IgG	Calculated MW: 510 aa, 59 kDa	
Immunogen Catalog Number: AG22802	Observed MW: 59 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : A549 cells, Jurkat cells
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 25724-1-AP (ZNF419 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.