

For Research Use Only

# GLI1 Polyclonal antibody

Catalog Number: 25733-1-AP

12 Publications



## Basic Information

Catalog Number:

25733-1-AP

Size:

150ul, Concentration: 260 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22962

GenBank Accession Number:

BC013000

GeneID (NCBI):

2735

UNIPROT ID:

P08151

Full Name:

GLI family zinc finger 1

Calculated MW:

1106 aa, 118 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

human, mouse

## Notable Publications

Author	Pubmed ID	Journal	Application
Yu Xu	36247302	Am J Transl Res	WB
Jilian R Melamed	29930746	Oncotarget	IF
Chia-Hao Huang	27377153	Lab Chip	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

