

For Research Use Only

OATP5A1/SLCO5A1 Polyclonal antibody, PBS Only

Catalog Number: 25744-1-PBS



Basic Information

Catalog Number:

25744-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22620

GenBank Accession Number:

BC137424

GeneID (NCBI):

81796

UNIPROT ID:

Q9H2Y9

Full Name:

solute carrier organic anion transporter family, member 5A1

Observed MW:

110-12 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25744-1-AP (SLCO5A1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25744-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 25744-1-AP (SLCO5A1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25744-1-PBS in a different storage buffer formulation.