

For Research Use Only

P2RX2 Polyclonal antibody

Catalog Number: 25790-1-AP



Basic Information

Catalog Number:

25790-1-AP

Size:

150ul , Concentration: 160 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22904

GenBank Accession Number:

BC109200

GeneID (NCBI):

22953

UNIPROT ID:

Q9UBL9

Full Name:

purinergic receptor P2X, ligand-gated ion channel, 2

Calculated MW:

497 aa, 55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

