

For Research Use Only

ZNF672 Polyclonal antibody, PBS Only

Catalog Number: 25794-1-PBS



Basic Information

Catalog Number: 25794-1-PBS	GenBank Accession Number: BC050348	Purification Method: Antigen affinity purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 79894	
Source: Rabbit	UNIPROT ID: Q499Z4	
Isotype: IgG	Full Name: zinc finger protein 672	
Immunogen Catalog Number: AG22886	Calculated MW: 452 aa, 50 kDa	
	Observed MW: 50-60 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

ZNF672 (Zinc Finger Protein 672) belongs to the Krueppel C2H2-type zinc-finger protein family, a large class of transcription factors characterized by tandem C2H2-type zinc finger domains that mediate sequence-specific DNA binding.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Rat liver tissue lysates were subjected to SDS PAGE followed by western blot with 25794-1-AP (ZNF672 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25794-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using ZNF672 antibody (25794-1-AP) at dilution of 1:400 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 25794-1-PBS in a different storage buffer formulation.