

For Research Use Only

SPATC1 Polyclonal antibody

Catalog Number: 25861-1-AP

2 Publications



Basic Information

Catalog Number:

25861-1-AP

Size:

150ul, Concentration: 550 ug/ml by Nanodrop and 307 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12157

GenBank Accession Number:

BC053547

GeneID (NCBI):

375686

UNIPROT ID:

Q76KD6

Full Name:

spermatogenesis and centriole associated 1

Calculated MW:

500 aa, 52 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB : human testis tissue,

Notable Publications

Author	Pubmed ID	Journal	Application
Brett Nixon	33451771	Mol Cell Proteomics	IF
Brett Nixon	30072580	Mol Cell Proteomics	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 25861-1-AP (SPATC1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.