

For Research Use Only

NUDT15 Polyclonal antibody

Catalog Number: 25914-1-AP



Basic Information

Catalog Number:

25914-1-AP

Size:

150ul, Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23242

GenBank Accession Number:

BC133015

GeneID (NCBI):

55270

UNIPROT ID:

Q9NV35

Full Name:

nudix (nucleoside diphosphate linked moiety X)-type motif 15

Observed MW:

19 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells, HeLa cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat lysates were subjected to SDS PAGE followed by western blot with 25914-1-AP (NUDT15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.