

For Research Use Only

HARBI1 Polyclonal antibody

Catalog Number: 25946-1-AP



Basic Information

Catalog Number:

25946-1-AP

Size:

150ul, Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22522

GenBank Accession Number:

BC036925

GeneID (NCBI):

283254

UNIPROT ID:

Q96MB7

Full Name:

harbinger transposase derived 1

Calculated MW:

349 aa, 39 kDa

Observed MW:

~39 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : RAW 264.7 cells, THP-1 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25946-1-AP (HARBI 1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.