

For Research Use Only

HARBI1 Polyclonal antibody, PBS Only

Catalog Number: 25946-1-PBS

Featured Product



Basic Information

Catalog Number:

25946-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22522

GenBank Accession Number:

BC036925

GeneID (NCBI):

283254

UNIPROT ID:

Q96MB7

Full Name:

harbinger transposase derived 1

Calculated MW:

349 aa, 39 kDa

Observed MW:

~39 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25946-1-AP (HARBI 1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25946-1-PBS in a different storage buffer formulation.



WB result of HARBI 1 antibody (25946-1-AP; 1:4000; incubated at room temperature for 1.5 hours) with sh-Control and sh-HARBI1 transfected Jurkat cells. This data was developed using the same antibody clone with 25946-1-PBS in a different storage buffer formulation.