

For Research Use Only

IL-6 Polyclonal antibody

Catalog Number: 26049-1-AP

3 Publications



Basic Information

Catalog Number:

26049-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23467

GenBank Accession Number:

NM_031168

GeneID (NCBI):

16193

Full Name:

interleukin 6

Calculated MW:

24 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, ELISA

Species Specificity:

human

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Zhigang Liu	36288667	Bioorg Chem	WB
Zhao Chen	36410585	Free Radic Biol Med	ELISA
Xin Yu	39913044	Biochem Genet	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

