

For Research Use Only

FLJ25006 Polyclonal antibody

Catalog Number: 26103-1-AP



Basic Information

Catalog Number: 26103-1-AP	GenBank Accession Number: BC045622	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 450 µg/ml by Nanodrop and 347 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 124923	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: uncharacterized serine/threonine-protein kinase SgK494	
Isotype: IgG	Calculated MW: 33 kDa	
Immunogen Catalog Number: AG22017	Observed MW: 33 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 26103-1-AP (FLJ25006 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.