

For Research Use Only

DMRTB1 Polyclonal antibody

Catalog Number: 26118-1-AP



Basic Information

Catalog Number: 26118-1-AP	GenBank Accession Number: BC029566	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 63948	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: DMRT-like family B with proline-rich C-terminal, 1	
Isotype: IgG	Observed MW: 40-45 kDa	
Immunogen Catalog Number: AG23639		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human testis tissue,
Species Specificity: human	

Background Information

DMRTB1, also named as Doublesex- and mab-3-related transcription factor B1, is a 342 amino acid protein, which contains 1 DM DNA-binding domain and belongs to the DMRT family. DMRTB1 is expressed in testis.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 26118-1-AP (DMRTB1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.