

For Research Use Only

ECHDC2 Polyclonal antibody

Catalog Number: 26126-1-AP



Basic Information

Catalog Number:

26126-1-AP

Size:

150ul, Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23510

GenBank Accession Number:

BC044574

GeneID (NCBI):

55268

UNIPROT ID:

Q86YB7

Full Name:

enoyl Coenzyme A hydratase domain containing 2

Observed MW:

29-31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : rat kidney tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat kidney tissue were subjected to SDS PAGE followed by western blot with 26126-1-AP (ECHDC2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.