

For Research Use Only

ECHDC2 Polyclonal antibody, PBS Only

Catalog Number: 26126-1-PBS



Basic Information

Catalog Number: 26126-1-PBS	GenBank Accession Number: BC044574	Purification Method: Antigen affinity purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 55268	
Source: Rabbit	UNIPROT ID: Q86YB7	
Isotype: IgG	Full Name: enoyl Coenzyme A hydratase domain containing 2	
Immunogen Catalog Number: AG23510	Observed MW: 29-31 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, rat

Background Information

ECHDC2 (Enoyl-CoA Hydratase Domain Containing 2) is a protein-coding gene that encodes a mitochondrial enzyme. ECHDC2 is involved in fatty acid β -oxidation and branched-chain amino acid catabolism.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat kidney tissue were subjected to SDS PAGE followed by western blot with 26126-1-AP (ECHDC2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26126-1-PBS in a different storage buffer formulation.