

For Research Use Only

C14orf28 Polyclonal antibody, PBS Only

Catalog Number: 26130-1-PBS



Basic Information

Catalog Number:

26130-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23441

GenBank Accession Number:

BC105030

GeneID (NCBI):

122525

UNIPROT ID:

Q4W4Y0

Full Name:

chromosome 14 open reading frame 28

Observed MW:

35-40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, IP, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 26130-1-AP (C14orf28 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26130-1-PBS in a different storage buffer formulation.



IP result of anti-C14orf28 (IP:26130-1-AP, 4ug; Detection:26130-1-AP 1:300) with HeLa cells lysate 1200ug. This data was developed using the same antibody clone with 26130-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using C14orf28 antibody (26130-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 26130-1-PBS in a different storage buffer formulation.