

For Research Use Only

SPRR2A Polyclonal antibody

Catalog Number: 26135-1-AP



Basic Information

Catalog Number: 26135-1-AP	GenBank Accession Number: BC096108	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 300 ug/ml by Nanodrop;	GeneID (NCBI): 6700	Recommended Dilutions: IHC: 1:50-1:500
Source: Rabbit	UNIPROT ID: P35326	
Isotype: IgG	Full Name: small proline-rich protein 2A	
Immunogen Catalog Number: AG23537		

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human colon cancer tissue,
Species Specificity: human	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 26135-1-AP (SPRR2A antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).