

For Research Use Only

# DUB3 Polyclonal antibody

Catalog Number: 26143-1-AP

Featured Product

8 Publications



## Basic Information

Catalog Number:

26143-1-AP

Size:

150ul, Concentration: 280 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23567

GenBank Accession Number:

BC157073

GeneID (NCBI):

377630

UNIPROT ID:

Q6R6M4

Full Name:

deubiquitinating enzyme 3

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC, IP

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB : SGC-7901 cells,

## Notable Publications

Author	Pubmed ID	Journal	Application
Haile Zhao	34658851	Front Pharmacol	WB,IP
Zhitao Jiang	35849971	Phytomedicine	WB
Yang Su	36040642	Breast Cancer Res Treat	WB,IP

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26143-1-AP (DUB3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.