

For Research Use Only

DUB3 Polyclonal antibody

Catalog Number: 26143-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number: 26143-1-AP	GenBank Accession Number: BC157073	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 280 µg/ml by Nanodrop;	GeneID (NCBI): 377630	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: deubiquitinating enzyme 3	
Isotype: IgG	Observed MW: 60 kDa	
Immunogen Catalog Number: AG23567		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : SGC-7901 cells,
Cited Applications: IP, WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Haile Zhao	34658851	Front Pharmacol	WB,IP
Zhitao Jiang	35849971	Phytomedicine	WB
Yang Su	36040642	Breast Cancer Res Treat	WB,IP

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26143-1-AP (DUB3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.