

For Research Use Only

# Cytokeratin 71 Polyclonal antibody, PBS Only

Catalog Number: 26149-1-PBS



## Basic Information

<b>Catalog Number:</b> 26149-1-PBS	<b>GenBank Accession Number:</b> BC103917	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 112802	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q3SY84	
<b>Isotype:</b> IgG	<b>Full Name:</b> keratin 71	
<b>Immunogen Catalog Number:</b> AG23536	<b>Observed MW:</b> 57 kDa	

## Applications

**Tested Applications:**  
WB, Indirect ELISA

**Species Specificity:**  
human, mouse

## Background Information

KRT71 belongs to a family of type II keratins that are specifically expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q13.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 26149-1-AP (KRT71 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26149-1-PBS in a different storage buffer formulation.