

For Research Use Only

PATZ1 Polyclonal antibody

Catalog Number: 26225-1-AP

1 Publications



Basic Information

Catalog Number:

26225-1-AP

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop and 527 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23912

GenBank Accession Number:

BC021091

GeneID (NCBI):

23598

UNIPROT ID:

Q9HBE1

Full Name:

POZ (BTB) and AT hook containing zinc finger 1

Observed MW:

70-74 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse spleen tissue, Raji cells

Notable Publications

Author	Pubmed ID	Journal	Application
Cuihong Zhao	29926841	Med Sci Monit	WB,IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 26225-1-AP (PATZ1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.