

For Research Use Only

PDZD11 Polyclonal antibody, PBS Only

Catalog Number: 26236-1-PBS



Basic Information

Catalog Number:

26236-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23975

GenBank Accession Number:

BC012996

GeneID (NCBI):

51248

UNIPROT ID:

Q5EBL8

Full Name:

PDZ domain containing 11

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

PDZ Domain Containing 11 (PDZD11), a small protein widely expressed in the body, is mainly comprised of a single PDZ domain. PDZD11 has often been reported to be associated with the adherens junction (AJ) of epithelial cells. PDZD11 binds nectins to the cadherin-and cytoskeleton-related protein complex by interacting with PLEKHA7 to stabilize the nexins at AJs and promotes efficient early assembly of junctions. (PMID: 38766598, PMID: 35714771)

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 26236-1-AP (PDZD11 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26236-1-PBS in a different storage buffer formulation.