

For Research Use Only

ZXDA Polyclonal antibody

Catalog Number: 26261-1-AP



Basic Information

Catalog Number: 26261-1-AP	GenBank Accession Number: BC059356	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 7789	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: zinc finger, X-linked, duplicated A	
Isotype: IgG	Observed MW: 85-90 kDa	
Immunogen Catalog Number: AG23411		

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : A431 cells,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 26261-1-AP (ZKDA Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.