

For Research Use Only

WDR59 Polyclonal antibody

Catalog Number: 26289-1-AP



Basic Information

Catalog Number:

26289-1-AP

Size:

150ul , Concentration: 750 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23495

GenBank Accession Number:

BC014887

GeneID (NCBI):

79726

UNIPROT ID:

Q6PJI9

Full Name:

WD repeat domain 59

Observed MW:

110 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human, mouse

Positive Controls:

IP : mouse brain tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-WDR59 (IP:26289-1-AP, 4ug;
Detection:26289-1-AP 1:800) with mouse brain
tissue lysate 1840 ug.