

For Research Use Only

CCDC135 Polyclonal antibody

Catalog Number: 26307-1-AP



Basic Information

Catalog Number:

26307-1-AP

Size:

150ul, Concentration: 680 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23635

GenBank Accession Number:

BC036667

GeneID (NCBI):

84229

UNIPROT ID:

Q8IY82

Full Name:

coiled-coil domain containing 135

Observed MW:

103 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human testis tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26307-1-AP (CCDC135 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.