

For Research Use Only

CCDC135 Polyclonal antibody, PBS Only

Catalog Number: 26307-1-PBS



Basic Information

Catalog Number: 26307-1-PBS	GenBank Accession Number: BC036667	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 84229	
Source: Rabbit	UNIPROT ID: Q8IY82	
Isotype: IgG	Full Name: coiled-coil domain containing 135	
Immunogen Catalog Number: AG23635	Observed MW: 103 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26307-1-AP (CCDC135 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26307-1-PBS in a different storage buffer formulation.