

For Research Use Only

RFX2 Polyclonal antibody

Catalog Number: 26313-1-AP



Basic Information

Catalog Number: 26313-1-AP	GenBank Accession Number: BC028579	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 550 µg/ml by Nanodrop;	GeneID (NCBI): 5990	Recommended Dilutions: IP 0.5-4.0 ug for IP and 1:200-1:1000 for WB
Source: Rabbit	Full Name: regulatory factor X, 2 (influences HLA class II expression)	
Isotype: IgG	Calculated MW: 723 aa, 80 kDa	
Immunogen Catalog Number: AG24737	Observed MW: 80 kDa	

Applications

Tested Applications: IP, ELISA	Positive Controls: IP : mouse testis tissue,
Species Specificity: mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-RFX2(IP:26313-1-AP, 4ug;
Detection:26313-1-AP 1:400) with mouse testis
tissue lysate 1120 ug.