

For Research Use Only

# RFX2 Polyclonal antibody

Catalog Number: 26313-1-AP



## Basic Information

### Catalog Number:

26313-1-AP

### Size:

150ul , Concentration: 550 ug/ml by Nanodrop;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG24737

### GenBank Accession Number:

BC028579

### GeneID (NCBI):

5990

### UNIPROT ID:

P48378

### Full Name:

regulatory factor X, 2 (influences HLA class II expression)

### Calculated MW:

723 aa, 80 kDa

### Observed MW:

80 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

## Applications

### Tested Applications:

IP, ELISA

### Species Specificity:

mouse

### Positive Controls:

IP : mouse testis tissue,

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



IP result of anti-RFX2 (IP:26313-1-AP, 4ug;  
Detection:26313-1-AP 1:400) with mouse testis  
tissue lysate 1120 ug.