

For Research Use Only

GTSF1 Polyclonal antibody

Catalog Number: 26340-1-AP



Basic Information

Catalog Number: 26340-1-AP	GenBank Accession Number: BC021179	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 121355	Recommended Dilutions: WB 1:5000-1:50000
Source: Rabbit	Full Name: gametocyte specific factor 1	
Isotype: IgG	Observed MW: 19-23 kDa	
Immunogen Catalog Number: AG23647		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : K-562 cells, mouse testis tissue, rat testis tissue
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26340-1-AP (GTSF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.