

For Research Use Only

NDUFAF5 Polyclonal antibody

Catalog Number: 26359-1-AP



Basic Information

Catalog Number: 26359-1-AP	GenBank Accession Number: BC005984	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 79133	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: chromosome 20 open reading frame 7	
Isotype: IgG	Observed MW: 29 kDa	
Immunogen Catalog Number: AG23822		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 26359-1-AP (NDUFAF5 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.