

For Research Use Only

ZNF16 Polyclonal antibody, PBS Only

Catalog Number: 26374-1-PBS



Basic Information

Catalog Number:

26374-1-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23835

GenBank Accession Number:

BC010996

GeneID (NCBI):

7564

UNIPROT ID:

P17020

Full Name:

zinc finger protein 16

Observed MW:

70-77 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Background Information

ZNF16 (also known as HZF1 and KOX9) is a multi-C2H2 zinc finger protein first identified via its expression in human T-cells and shown to have a role in blood cell differentiation. ZNF16 was later shown to be ubiquitously expressed in a variety of fetal and adult tissues, suggesting a broader function. (PMID: 41279147, PMID: 41666695)

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26374-1-AP (ZNF16 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26374-1-PBS in a different storage buffer formulation.