

For Research Use Only

FOPNL Polyclonal antibody

Catalog Number: 26409-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number: 26409-1-AP	GenBank Accession Number: BC022321	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 123811	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: chromosome 16 open reading frame 63	
Isotype: IgG	Observed MW: 22 kDa	
Immunogen Catalog Number: AG23991		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells,
Cited Applications: IF, IHC	
Species Specificity: human	
Cited Species: zebrafish	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shanshan Xie	30475641	FASEB J	IF,IHC

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 26409-1-AP (FOPNL Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.