

For Research Use Only

Macoilin Polyclonal antibody

Catalog Number: 26488-1-AP



Basic Information

Catalog Number:

26488-1-AP

Size:

150ul, Concentration: 850 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24428

GenBank Accession Number:

BC032427

GeneID (NCBI):

55219

UNIPROT ID:

Q8N5G2

Full Name:

transmembrane protein 57

Calculated MW:

664 aa, 76 kDa

Observed MW:

76 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP: Neuro-2a cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-Macoilin (IP:26488-1-AP, 4ug;
Detection:26488-1-AP 1:2000) with Neuro-2a cells
lysate 1440 ug.