

For Research Use Only

Macoilin Polyclonal antibody, PBS Only

Catalog Number: 26488-1-PBS



Basic Information

Catalog Number:

26488-1-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24428

GenBank Accession Number:

BC032427

GeneID (NCBI):

55219

UNIPROT ID:

Q8N5G2

Full Name:

transmembrane protein 57

Calculated MW:

664 aa, 76 kDa

Observed MW:

76 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-Maccoilin (IP:26488-1-AP, 4ug; Detection:26488-1-AP 1:2000) with Neuro-2a cells lysate 1440 ug. This data was developed using the same antibody clone with 26488-1-PBS in a different storage buffer formulation.