

For Research Use Only

FAM161B Polyclonal antibody

Catalog Number: 26492-1-AP



Basic Information

Catalog Number: 26492-1-AP	GenBank Accession Number: BC053909	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 450 ug/ml by Nanodrop;	GeneID (NCBI): 145483	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q96MY7	
Isotype: IgG	Full Name: family with sequence similarity 161, member B	
Immunogen Catalog Number: AG24566	Observed MW: 72-74 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : THP-1 cells, U-251 cells
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26492-1-AP (FAM161B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.