

For Research Use Only

FAM161B Polyclonal antibody, PBS Only

Catalog Number: 26492-1-PBS



Basic Information

Catalog Number: 26492-1-PBS	GenBank Accession Number: BC053909	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 145483	
Source: Rabbit	UNIPROT ID: Q96MY7	
Isotype: IgG	Full Name: family with sequence similarity 161, member B	
Immunogen Catalog Number: AG24566	Observed MW: 72-74 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26492-1-AP (FAM161B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26492-1-PBS in a different storage buffer formulation.